

ANNEX

DETERMINATION OF DICHLOROMETHANE AND 1,1,1-TRICHLOROETHANE QUANTITATIVE DETERMINATION OF TOSYLCHLORAMIDE SODIUM (INN) (CHLORAMINE-T)

6. PROCEDURE

6.2. *Extraction*

- 6.2.1. Take 30 ml of the filtrate and extract three times with 15 ml of diethyl ether (4.4). If necessary dry the ethereal phases and collect them in a 50 ml graduated flask. Make up with diethyl ether (4.4).
- 6.2.2. Take 25 ml of the dried ethereal extract and evaporate to dryness in a nitrogen stream. Redissolve the residue with 1 ml of ethanol (4.5).